

New diagnostic tools for improved diagnostics of grapevine phytoplasmas


Mehle Nataša, Ravnikar Maja, Kogovšek Polona, Jakomin Tjaša, Pugelj Anja, Dermastia Marina

natasa.mehle@nib.si

Phytoplasma

- cell wall-less Gram positive bacteria
- class Mollicutes
- cell and genome size are the smallest among bacteria
- obligate intracellular parasites
- Transmitted:
 - phloem-feeding leafhoppers, planthoppers and psyllids
 - dodder, micropropagation, grafting and cutting
- >1000 diseases

phytoplasmas in
phloem sieve
element



Grapevine yellows

- caused by different phytoplasmas (different vectors)
- indistinguishable by symptoms




Slovenia:

- '*Ca. P. solani*' -> **BN**

- Flavescence dorée phytoplasma -> **FD**

- '*Ca. P. asteris*' -> **AY**




2005-2015 (2234 samples)

Limitations of phytoplasma detection

- the smallest by size and genome
- routinely uncultivable – ~~traditional diagnostic methods suitable for bacteria~~
- uneven distribution in the phloem (vascular tissue in stem, leaves, roots)
- low concentration
- variations in titer according to the season/plant organ



FD is listed in the EU2000/29 Council Directive on Harmful organisms and the A2 quarantine list of pests of EPPO: **the destruction of diseased stocks, plants showing symptoms and surrounding plants is mandatory.**

Example: FD – Izola



Reliable, sensitive and fast diagnostic procedure is needed!

Diagnostic procedure



Phytoplasma

Methods and Protocols

Series: » Methods in Molecular Biology, Vol. 938

Dickinson, Matt; Hodggets, Jennifer (Eds.)

2013, 2013, XIII, 421 p. 68 illus., 45 in color.


A product of Humana Press

- | | | |
|----|--|-----|
| 12 | Automated DNA Extraction for Large Numbers of Plant Samples | 139 |
| | <i>Nataša Mehle, Petra Nikolić, Matevž Rupar, Jana Boben, Maja Ravnikar, and Marina Dermastia</i> | |
| 22 | A Real-Time PCR Detection System for the Bois Noir and Flavescence Dorée Phytoplasmas and Quantification of the Target DNA | 253 |
| | <i>Nataša Mehle, Nina Prezelj, Matjaž Hren, Jana Boben, Kristina Gruden, Maja Ravnikar, and Marina Dermastia</i> | |

The validation data about this method is available at EPPO website: <http://dc.eppo.int/validationlist.php>



Diagnostic procedure







P	F	K	qPCR	D
	P	F		

+ less contamination, higher sensitivity

P	N2	CTAB	AGE	3x PCR	AGE	D	2x nPCR	AGE	D	nPCR	AGE	RLFP	PAGE	D
---	----	------	-----	--------	-----	---	---------	-----	---	------	-----	------	------	---



1st day	2nd day	3rd day	4th day	5th day
---------	---------	---------	---------	---------




LAMP : Loop mediated isothermal AMPLification

Fast multiplication of DNA targets in one tube
 at isothermal conditions (60-65°C)
 using a set of 4 or 6 primers

-Relatively simple application
 -No expensive equipment needed




Starting structure producing step




Cycling amplification step

Detection of LAMP products

- LAMP product on gel




- Turbidity




- LFD



- Fluorescence





- Array tubes



- Real time

- Intercalating dye (!)
- Fluorescent probes




Simultaneously heater and fluorimeter (e.g., GenieII/III, SmartCycler)

Real time detection of LAMP products


1) READING RESULTS

pos: rise of fluorescence
neg: no rise of fluorescence



2) CONFIRMATION OF RESULTS

Melting temperature of the final product is pathogen specific



Legend:

- positive control (amplification-> rise of fluorescence)
- sample (comparable to positive control)
- negative control (no amplification-> no fluorescence)

LAMP detection of phytoplasmas FD and BN

- FDp:

Plant Pathology (2015) 64, 286–296

Doi: 10.1111/ppa.12266



LAMP assay and rapid sample preparation method for on-site detection of flavescence dorée phytoplasma in grapevine

P. Kogovšek^{ab*}, J. Hodgetts^c, J. Hall^c, N. Prezelj^a, P. Nikolić^a, N. Mehle^a, R. Lenarčič^a, A. Rotter^a, M. Dickinson^d, N. Boonham^c, M. Dermastia^a and M. Ravnikar^a

^aDepartment of Biotechnology and Systems Biology, National Institute of Biology; ^bDepartment of Biology, Biotechnical Faculty, University of Ljubljana, Večna pot 111, 1000 Ljubljana, Slovenia; ^cThe Food and Environment Research Agency, Sand Hutton, York YO41 1LZ; and ^dSchool of Biosciences, University of Nottingham, Sutton Bonington Campus, Loughborough LE12 5RD, UK



-BNp:

Euphresco GRAFDEPI 2

On-site application of the FDp and BNp testing



Sampling



Homogenisation



No need for DNA extraction and purification (not sensitive to inhibitors)




Mix the sample with LAMP reagents




Amplification & Detection

Comparison of time needed for FDp and BNP detection with different methods



LAMP – validation (FD example)




Dilution	FDp DNA copy no.	Extraction of DNA with KingFisher + qPCR (Cq)	LAMP (Tp)
3x	243-729	+ (27.9)	+ (21.1)
9x	81-243	+ (29.5)	+ (27.3)
27x	27-81	+ (31.4)	+ (25.0)
81x	9-27	+ (32.9)	+ (19.1)
243x	3-9	+ (34.4)	-
729x	1-3	+ (34.8)	-
2187x	0	-	-

LAMP is 9x less sensitive than qPCR (analytical sensitivity)

(Kogovšek et al., 2015)

LAMP – validation (FD example)



Diagnostic sensitivity:


**-Testing of extracted DNA (52
FDp infected samples): 100%**

**-Testing of crude homogenates
(27 FDp infected samples): 100%**

LAMP – validation

(FD & BN – testing of crude homogenates)

2015:
286 official
grapevine samples



		Cq < 32	Cq > 32		
BN pos	212	193	19	12	5
				2	
FD pos		Cq < 31	Cq > 31		
	28	24	4	3	1

LAMP – test performance study

(FD & BN – testing of extracted DNA)

- Euphresco project GRAFDEPI 2
- Participants: 10 laboratories (from the research and plant protection area from Europe and Australia)
- Additionally, LAMP FDp assay (Kogovšek et al., 2015) was compared with a Qualiplante/Hyris isothermal amplification assay for FD (ISOA FD Qualiplante) by 3 laboratories
- 18 DNA samples were subject of this TPS

LAMP – test performance study

(FD & BN – testing of extracted DNA)

	Assay		
	LAMP BN	LAMP FD (Kogovšek et al., 2015)	ISOA FD Qualiplante
No. of labs taking in account for the evaluation	10	10	3
No. of results	180	179	54
N ⁺	50	49	15
PA	49	49	15
ND	1	0	0
Undetermined (sus) of N ⁺	0	0	0
N ⁻	130	130	39
NA	130	127	38
PD	0	0	0
Undetermined (sus) of N ⁻	0	3	1
Accuracy	99,4%	98,3%	98,1%
Rate of true positives	98,0%	100%	100%
Rate of true negatives	100%	97,7%	97,4%

The validation data for LAMP FD (Kogovšek et al., 2015) available at EPPO website:
<http://dc.eppo.int/validationlist.php>

Quantification

- monitor phytoplasma kinetics (progress of an infection, and variations of the phytoplasma titer through the season and in different plant tissues)



Plant Pathology (2013)

Doi: 10.1111/j.1365-3059.2012.02693.x


Spatiotemporal distribution of flavescence dorée phytoplasma in grapevine

N. Prezelj, P. Nikolić, K. Gruden, M. Ravnikar and M. Dermastia*

- screening plants for resistance against phytoplasma
- estimate the number of copies carried by the vectors

Quantification

- Real time PCR:
quantification
against reference
material (standard
curve):



No certified phytoplasma reference material (dilutions of a sample containing the target DNA sequence or a sample with known copy numbers of plasmids)

Quantification


- **Digital PCR**

- absolute quantification of target sequences without relying on the use of standard curves


- **droplet digital PCR (ddPCR):**



droplet generation



amplification (PCR)



reading

Analysis of FDp with ddPCR

- Transfer from qPCR to ddPCR

Plant Pathology (2007) 56, 785–796

Doi: 10.1111/j.1365-3059.2007.01688.x

Real-time PCR detection systems for *Flavescence dorée* and *Bois noir* phytoplasmas in grapevine: comparison with conventional PCR detection and application in diagnostics

M. Hren^{a*}, J. Boben^a, A. Rotter^a, P. Kralj^b, K. Gruden^a and M. Ravnikar^a

- same primers and probes
- change in mastermix

Phytopathogenic Mollicutes
Vol. 4(1), June 2014, 9-15

doi : 10.5958/2249-4677.2014.00576.3

IndianJournals.com
A product of Diva Enterprises Pvt. Ltd.



Research Article

Quantitative analysis of “flavescence doreé” phytoplasma with droplet digital PCR



Nataša Mehle, Tanja Dreo and Maja Ravnikar

Analysis of FDp with ddPCR

Sensitivity: similar as with qPCR


qPCR

ddPCR




Analysis of FDp with ddPCR

Repeatability of **ddPCR** and **qPCR**:



higher precision and repeatability of ddPCR for quantification of FDp at the low concentrations



Conclusions – phytoplasma detection

- **Diagnostic procedure:**

simple&quick homogenisation step + DNA extraction based on the binding of DNA to magnetic beads + real-time PCR


- **LAMP assay for FDp and BNp:**

- Application in laboratories (high through-put) or without expensive equipment on-site
- LAMP is less prone to inhibition therefore just homogenization of samples without NA extraction is sufficient

Conclusions – phytoplasma detection

- **ddPCR for FDp:**

- Absolutely quantify phytoplasma without the need of any calibrant (calibration curves for quantification of FDp are not needed)
- Quantification and quality control of DNA based on in-house reference materials typically used in diagnostics and metrological laboratories



Descriptive assessment of uncertainties of qualitative real-time PCR for detection of plant pathogens and quality performance monitoring

N. Mehle¹, T. Dreo¹, C. Jeffries² and M. Ravnikar¹

¹National Institute of Biology, Večna pot 111, SI-1000, Ljubljana, Slovenia; e-mail: Natasa.mehle@nib.si

²Science and Advice for Scottish Agriculture, Roddinglaw Road, EH12 9FJ, Edinburgh, UK

Acknowledgement

- FP7 project Vitisens
- Slovenian Research Agency (contract no. L4-5525)
- Slovenian phytosanitary administration

Acknowledgement

- Euphresco project GRAFDEPI 2

Acronim	Contact	Laboratory
AU-AgriBio	Fiona Constable	AgriBio, Australia
BE-ILVO	Kris De Jonghe	Institute for Agricultural and Fisheries Research (ILVO), Plant Sciences – Crop Protection
FR-ANSES	Marianne Loiseau	Laboratoire de la santé des végétaux / Plant Health Laboratory
IT-CNR	Flavio Veratti, Cristina Marzachi	Istituto per la Protezione Sostenibile delle Piante, CNR
IT-CRA-PAV	Luca Ferretti	CRA-PAV, Centro di Ricerca per la Patologia Vegetale
IT-CRA-VIT	Elisa Angelini Luisa Filippin	C.R.A. - VIT Centro di Ricerca per la Viticoltura
PT-INIAV	Esmeraldina Sousa, Eugénia Andrade	Instituto Nacional de Investigação Agrária e Veterinária, I.P.; Unidade Estratégica de Investigação e Serviços de Sistemas Agrários e Florestais e Sanidade Vegetal
SI-NIB	Marina Dermastia, Nataša Mehle, Polona Kogovšek	National Institute of Biology
UK-FERA	Jennifer Hodgetts	Fera Science Ltd. (Fera)
UK-Nottingham	Matt Dickinson	School of Biosciences, University of Nottingham, Sutton Bonington Campus

